

## Product Datasheet

**Product code:** **K4MPTSP**

**Product type:** Stainless steel probes

**Made in:** Pakistan

**Description:** The K4MPTSS Stainless Steel spatula kit contains the MPTSS1, MPTSS2, MPTSS3 and MPTSS4 models.

**Body material:** Stainless Steel (S)

**Picture:**



<b>Product dimensions:</b>	OAL	190.0 mm – 7.480 inch
	Width	70.0 mm – 2.756 inch
	Height	7.0 mm – 0.276 inch