


Cobalt Alloy KO Biology Tweezers

We added to our Tweezer portfolio 4 new high precision tweezers, made of Cobalt Alloy, meticulously crafted for demanding applications.



COBALT ALLOY MODELS

3.KO

4.KO

5.KO

7.KO

Precision for Critical Sectors

Whether it's aerospace, medical device instruments, or electronics, our tweezers are tailored for industries where there's no room for error.

Unparalleled Durability

Made with Cobalt Alloy (KO), these tweezers offer unmatched fatigue resistance and elasticity. They're built to last, ensuring consistent precision throughout their lifespan.

Popular Styles, Now in KO

Each of our new models - style 3, 4, 5, and 7 - boasts a distinct tip geometry, ensuring professionals have the exact tool tailored for their specific tasks.

Cobalt Alloy KO Biology Tweezers



3.KO

Tips: straight, ultra fine, pointed, superior finish.
OAL: 120mm - 4.72"



4.KO

Tips: straight, ultra fine, sharp, superior finish.
OAL: 110mm - 4.33"



5.KO

Tips: straight, ultra fine, sharp, superior finish.
OAL: 110mm - 4.33"



7.KO

Tips: curved, ultra fine, superior finish.
OAL: 115mm - 4.53"