


# MPT Stainless Steel probes dimensions



**MPTSP1** - OAL: 155 mm

	WIDTH (mm)	THICKNESS (mm)
TIP	0.25	0.25

**MPTSP2** - OAL: 155 mm

	WIDTH (mm)	THICKNESS (mm)
TIP	0.20	0.20

**MPTSP3** - OAL: 155 mm

	WIDTH (mm)	THICKNESS (mm)
TIP	0.50	0.40

**MPTSP4** - OAL: 155 mm

	WIDTH (mm)	THICKNESS (mm)
TIP	0.50	0.40

**MPTSP5** - OAL: 155 mm

	WIDTH (mm)	THICKNESS (mm)
TIP	0.50	0.60

**MPTSP6** - OAL: 155 mm

	WIDTH (mm)	THICKNESS (mm)
TIP	2.20	0.70