

MPT Stainless Steel spatulas dimensions



MPTSS1

OAL: 180 mm

TIPS	WIDTH (mm)	THICKNESS (mm)
A	0.60	0.50
B	5.30	0.80

MPTSS2

OAL: 175 mm

TIPS	WIDTH (mm)	THICKNESS (mm)
A	13.25	0.50
B	6.00	1.20

MPTSS3

OAL: 165 mm

TIPS	WIDTH (mm)	THICKNESS (mm)
A	8.00	0.50
B	8.70	0.50

MPTSS4

OAL: 180 mm

TIPS	WIDTH (mm)	THICKNESS (mm)
A	10.00	0.30
B	0.60	0.20