

PHOTONIC

Photonic optics develops high quality LED illumination systems for many different applications, including microscopy. For over 30 years Photonic is providing modular systems at a good price/performance ratio, making it the ideal partner for Dino-Lite digital microscopes. Especially the Dino-Lite Long Working Distance and Extra Long Working Distance models can benefit from the flexible illumination systems of Photonic. It complements the integrated LED's of the Dino-Lite microscopes perfectly.

Desk set

The Photonic Desk-set is composed of the LED Hi-Power Spot (2 pieces), 2 spot holders (flexible arms) and 2 base plates.

The Hi-Power Spot must be attached to the included Spot holder. This Spot holder is then mounted onto the round Base plate. With this set you can create two individual Hi-power spots with flexible arms. As it is mounted on the baseplate, it can be placed anywhere around the sample creating total flexibility.

The brightness can be controlled by the compact control unit. Furthermore you can switch on/off the 2-arms individually.



Desk set

- High power LED spots, 2x
 - Homogeneous illumination
- 2-arms
 - Flexible arms
- 2 Base plates
 - Stable mounting on the base plate
- Control unit:
 - Brightness
 - on/off for each arm individually
- Easily placed around the sample:
- Works with almost all Dino-Lite digital microscope models



	PH-Desk set	
	Hi-Power spot arm, 2 pieces	Base plate, 2 pieces
ESD-Safe:	Yes with 4mm banana socket, plug not included	-
Dimensions:	55cm long, each arm	10cm diameter
Weight:	125g each	50g each
Compatible with:	All Dino-Lite models, except high magnification models	
Package contents:	Hi-Power Spot arm (2x), Spot holders (2x), Base plate (2x), Controller unit, power supply with EU & UK plug	
Warranty:	2 years of European warranty	

** All other Photonic items are available through Dino-Lite Europe / IDCP BV.*