

## PHOTONIC

Photonic optics develops high quality LED illumination systems for many different applications, including microscopy. For over 30 years Photonic is providing modular systems at a good price/performance ratio, making it the ideal partner for Dino-Lite digital microscopes. Especially the Dino-Lite Long Working Distance and Extra Long Working Distance models can benefit from the flexible illumination systems of Photonic. It complements the integrated LED's of the Dino-Lite microscopes perfectly.

### HPS set

The Photonic HPS-set is composed of the LED **Hi-Power Spot 2-arm** and the MS35/36B(E) adapters.

This spot is used for direct attachment to microscope stands, like the Dino-Lite MS35/36B(E) stands. As it is mounted on the stand, it is perfect for use in small spaces.

The brightness can be controlled by the compact control unit. Furthermore you can switch on/off the 2-arms individually.

The flexible arms can be easily adjusted to provide the best angle of illumination to your sample.



## HPS set

- High power LED spots
  - Life time of about 50.000 hours
  - No noise, no vibrations
- 2-arms
  - Flexible arms
  - Homogeneous illumination
- MS35/36B stand adapter
  - Stable mounting on the pillar
- Control unit:
  - Brightness
  - on/off for each arm individually
- Compact and space saving
- Works with almost all Dino-Lite digital microscope models



	PH-HPS set	
	Hi-Power spot 2-arm	Adapters MS35/36B(E)
<b>ESD-Safe:</b>	Yes with 4mm banana socket, plug not included	-
<b>Dimensions:</b>	55cm long, each arm	-
<b>Weight:</b>	650g	50g
<b>Compatible with:</b>	All Dino-Lite models with MS35/36B(E) stand, except high magnification models	
<b>Package contents:</b>	Hi-Power Spot 2-arms, Controller unit, power supply with EU & UK plug, MS35/36B(E) adapter	
<b>Warranty:</b>	2 years of European warranty	

\* All other Photonic items are available through Dino-Lite Europe / IDCP BV.