

Product Datasheet

Product code: **1.SA.GP.0**

Product type: Gold plated tweezers

Made in: Switzerland

Description: Type 1 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The coating of 2-microns thickness, made of pure 24-carat gold plating, provides excellent scratch resistance and makes them particularly suitable for handling nickel grids in immuno-gold work. The tweezers exhibit remarkable durability under high-volume autoclaving conditions, and the gold plating enhances their resistance to chemical corrosion, oxidation, and electrical conductivity.

Tip description: straight, fine, pointed, superior finish

Coating/finishing: Gold plated (GP)

Body material: Anti-Magnetic Anti-Acid Stainless Steel (SA)


Picture:



Tweezer dimensions:

| | |
|--------|-----------------------|
| OAL | 120.0 mm – 4.724 inch |
| Width | 10.0 mm – 0.394 inch |
| Height | 11.0 mm – 0.433 inch |

Tip details:



Tip dimensions:

| | |
|---------------|----------------------|
| Tip width | 0.15 mm – 0.006 inch |
| Tip thickness | 0.10 mm – 0.004 inch |