

Product Datasheet

Product code: **2A.SA.GP.0**

Product type: Gold plated tweezers

Made in: Switzerland

Description: Type 2A high precision tweezers are suitable to delicate applications which require holding a flat component firmly. The coating of 2-microns thickness, made of pure 24-carat gold plating, provides excellent scratch resistance and makes them particularly suitable for handling nickel grids in immuno-gold work. The tweezers exhibit remarkable durability under high-volume autoclaving conditions, and the gold plating enhances their resistance to chemical corrosion, oxidation, and electrical conductivity.

Tip description: straight, fine, round, flat, superior finish

Coating/finishing: Gold plated (GP)

Body material: Anti-Magnetic Anti-Acid Stainless Steel (SA)


Picture:




Tweezer dimensions:

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	11.0 mm – 0.433 inch

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

Tip width	1.90 mm – 0.075 inch
Tip thickness	0.11 mm – 0.004 inch