

## Product Datasheet

**Product code:** **2A.SA.T.0**

**Product type:** PTFE coated tweezers

**Made in:** Switzerland

**Description:** Type 2A high precision tweezers are suitable to delicate applications which require holding a flat component firmly. The PTFE (a material equivalent to Teflon®) coating provides the tweezers with chemical resistance properties suitable to acid and basic environments and for the cryogenic industry. The PTFE coating is used where low friction is required.

**Tip description:** straight, very fine, round, flat, superior finish

**Coating/finishing:** PTFE (T)

**Body material:** Anti-Magnetic Anti-Acid Stainless Steel (SA)


**Picture:**



**Tweezer dimensions:**

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	11.0 mm – 0.433 inch

**Tip details:**



**Tip dimensions:**

Tip width	1.90 mm – 0.075 inch
Tip thickness	0.11 mm – 0.004 inch