

Product Datasheet

Product code: **2WFSVR.SA.1**

Product type: Plastic Replaceable Tip Tweezers

Made in: Switzerland

Description: Type 2WFSV Plastic replaceable tip tweezers are suitable to delicate applications which require holding a 2-inch flat component firmly. The Carbon Fiber reinforced ESD PVDF Plastic provides the tweezers with ESD and exceptional non-scratching properties combined with excellent chemical resistance.

Tip description: straight, wafer

Tip material: PVDF (SV)

Coating/finishing: Sandblasted

Body material: Anti-Magnetic Anti-Acid Stainless Steel (SA)


Picture:




Tweezer dimensions:

OAL	130.0 mm – 5.118 inch
Width	12.0 mm – 0.472 inch
Height	18.0 mm – 0.709 inch

Tip details:



Dimensions inch



Dimensions mm

Tip dimensions:

Tip width	6.60 mm – 0.260 inch
B	10.50 mm – 0.413 inch
C	2.00 mm – 0.079 inch
D	3.50 mm – 0.138 inch
E	4.50 mm – 0.177 inch