

Product Datasheet

Product code: 3.KO.0

Product type: Biology Tweezers

Made in: Switzerland


Description: Type 3 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The Cobalt Alloy has a high fatigue resistance and maintain its elasticity throughout its lifespan, resulting in exceptional durability that is unmatched by other materials.

Tip description: straight, ultra fine, pointed, superior finish

Coating/finishing: Polished

Body material: Cobalt Alloy (KO)


Picture:



Tweezer dimensions:

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	9.0 mm – 0.354 inch

Tip details:



Tip dimensions:

Tip width	0.05 mm – 0.002 inch
Tip thickness	0.03 mm – 0.001 inch