

Product Datasheet

Product code: **315S-SC-90.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 315S-SC-90 are High precision scissors with long blades and long handles. The micro-serrated blade allows non-slip cutting. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, thick, round, micro-serrated

Coating/finishing: Mirror Polished

Handle coating: Black coating (BK)

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	230.0 mm – 9.055 inch
Width	60.0 mm – 2.362 inch
Height	7.5 mm – 0.295 inch

Tip dimensions:

Tip width	4.00 mm – 0.157 inch
Tip thickness	2.00 mm – 0.079 inch