

Product Datasheet

Product code: 330BMS.9

Product type: Surgical Scissors

Made in: Germany

Description: The 330B are Micro scissors with fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The curved blades prevent the tips from getting stuck in the material. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: curved, fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:	OAL	110.0 mm – 4.331 inch
	Width	24.0 mm – 0.945 inch
	Height	7.5 mm – 0.295 inch

Tip dimensions:	Tip width	0.30 mm – 0.012 inch
	Tip thickness	0.50 mm – 0.020 inch