

Product Datasheet

Product code: **331MS.9**

Product type: Surgical Scissors

Made in: Germany


Description: The 331 are Micro scissors with fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	110.0 mm – 4.331 inch
Width	24.0 mm – 0.945 inch
Height	7.5 mm – 0.295 inch

Tip dimensions:

Tip width	0.30 mm – 0.012 inch
Tip thickness	0.50 mm – 0.020 inch