

Product Datasheet

Product code: **335S-SC-67.P**

Product type: Surgical Scissors

Made in: PL

Description: The 335S-SC-67 are High precision scissors with long blades and long handles. The serrated edge provides the tips with extra grip, also preventing any possible slips. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, thick, round, polished

Coating/finishing: Sandblasted

Handle coating: Gold plated (GP)

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	170.0 mm – 6.693 inch
Width	60.0 mm – 2.362 inch
Height	6.5 mm – 0.256 inch

Tip dimensions:

Tip width	3.50 mm – 0.138 inch
Tip thickness	1.40 mm – 0.055 inch