

Product Datasheet

Product code: **343BMS.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 343BMS are high quality Micro scissors with very fine and accurate blades, ideal for precise and miniature work in confined spaces. The curved blades prevent the tips from getting stuck in the material. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The high-quality Stainless-steel material is suitable for cleaning. The blades are hardened.

Tip description: curved, very fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	105.0 mm – 4.134 inch
Width	25.0 mm – 0.984 inch
Height	7.0 mm – 0.276 inch

Tip dimensions:

Tip width	0.40 mm – 0.016 inch
Tip thickness	0.40 mm – 0.016 inch