

Product Datasheet

Product code: 345-55.B

Product type: Surgical Scissors

Made in: Pakistan

Description: The 345-55 are precision scissors with long blades and long handles. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: curved, fine, pointed

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:	OAL	145.0 mm – 5.709 inch
	Width	52.0 mm – 2.047 inch
	Height	13.0 mm – 0.512 inch

Tip dimensions:	Tip width	0.30 mm – 0.012 inch
	Tip thickness	1.00 mm – 0.039 inch