

Product Datasheet

Product code: **345BMS.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 345B are Micro scissors with very fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The curved blades prevent the tips from getting stuck in the material. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: curved, very fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	100.0 mm – 3.937 inch
Width	23.0 mm – 0.906 inch
Height	6.5 mm – 0.256 inch

Tip dimensions:

Tip width	0.30 mm – 0.012 inch
Tip thickness	0.50 mm – 0.020 inch