

Product Datasheet

Product code: **347MS.9**

Product type: Surgical Scissors

Made in: Germany


Description: The 347 are Micro scissors with very fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, very fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	100.0 mm – 3.937 inch
Width	24.0 mm – 0.945 inch
Height	7.0 mm – 0.276 inch

Tip dimensions:

Tip width	2.00 mm – 0.079 inch
Tip thickness	0.30 mm – 0.012 inch