

Product Datasheet

Product code: **351-35.B**

Product type: Surgical Scissors

Made in: Pakistan

Description: The 351 are precision scissors with extra fine blades, ideal for miniature work. The curved blades prevent the tips from getting stuck in the material. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: curved, very fine, sharp

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	90.0 mm – 3.543 inch
Width	45.0 mm – 1.772 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	0.30 mm – 0.012 inch
Tip thickness	0.60 mm – 0.024 inch