

## Product Datasheet

**Product code:** **351-45.P**

**Product type:** Surgical Scissors

**Made in:** PL

**Description:** The 351 are high precision scissors with extra fine blades, ideal for miniature work. The curved blades prevent the tips from getting stuck in the material. The Inox material is suitable for cleaning. The blades are hardened.

**Tip description:** curved, very fine, pointed

**Coating/finishing:** Sandblasted

**Body material:** Stainless Steel (Hardened AISI 420)

**Picture:**



**Product dimensions:**

OAL	110.0 mm – 4.331 inch
Width	45.0 mm – 1.772 inch
Height	3.0 mm – 0.118 inch

**Tip dimensions:**

Tip width	0.20 mm – 0.008 inch
Tip thickness	0.35 mm – 0.014 inch