

Product Datasheet

Product code: **351S-SCTC-45-2.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 351S-SCTC are High precision scissors with fine blades, ideal for miniature work. The Tungsten Carbide insert blades are designed to improve cutting performance and extend blade life. The serrated edge provides the tips with extra grip, also preventing any possible slips. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, micro-serrated

Tip material: Tungsten carbide (TX)

Coating/finishing: Polished

Handle coating: Black and Gold Plated handles (BK + GP)

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	115.0 mm – 4.528 inch
Width	48.0 mm – 1.890 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	0.50 mm – 0.020 inch
Tip thickness	1.20 mm – 0.047 inch