

Product Datasheet

Product code: 352S-35.B

Product type: Surgical Scissors

Made in: Pakistan

Description: The 352S are precision scissors with extra fine blades, ideal for miniature work. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, very fine, sharp

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	90.0 mm – 3.543 inch
Width	45.0 mm – 1.772 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	0.30 mm – 0.012 inch
Tip thickness	0.60 mm – 0.024 inch