

Product Datasheet

Product code: **352S-R40.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 352S-R are high precision scissors with extra fine blades, ideal for miniature work. Ribbon Style Finger Ring Handle enhances the lightness of the tool. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, very fine, sharp

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	105.0 mm – 4.134 inch
Width	52.0 mm – 2.047 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	0.20 mm – 0.008 inch
Tip thickness	0.70 mm – 0.028 inch