

Product Datasheet

Product code:	353S-47.9
Product type:	Surgical Scissors
Made in:	Germany
Description:	The 353S are High precision scissors with extra fine blades, ideal for miniature work. The Inox material is suitable for cleaning. The blades are hardened.
Tip description:	straight, fine, pointed
Coating/finishing:	Polished
Body material:	Stainless Steel (Hardened AISI 420)
Picture:	



Product dimensions:	OAL	120.0 mm – 4.724 inch
	Width	45.0 mm – 1.772 inch
	Height	5.0 mm – 0.197 inch
Tip dimensions:	Tip width	0.30 mm – 0.012 inch
	Tip thickness	1.00 mm – 0.039 inch