

Product Datasheet

Product code: 356S-40.B

Product type: Surgical Scissors

Made in: Pakistan

Description: The 356S are precision scissors with extra fine blades, ideal for miniature work. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, fine, beveled

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:	OAL	105.0 mm – 4.134 inch
	Width	46.0 mm – 1.811 inch
	Height	4.0 mm – 0.157 inch

Tip dimensions:	Tip width	1.50 mm – 0.059 inch
	Tip thickness	1.50 mm – 0.059 inch