

Product Datasheet

Product code: **358S-SU-35.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 358S-SU are high precision scissors with fine blades, ideal for specific work. The specific hook shape of the blade tip is suitable for removing suture threads. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, thick, round

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	90.0 mm – 3.543 inch
Width	42.5 mm – 1.673 inch
Height	3.5 mm – 0.138 inch

Tip dimensions:

Tip width	2.30 mm – 0.091 inch
Tip thickness	1.50 mm – 0.059 inch