

Product Datasheet

Product code: 366MS.B

Product type: Micro Scissors

Made in: Pakistan


Description: The 346AMS are Micro scissors with very fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The Inox material is suitable for cleaning. The blades are hardened to 55HRC hardness.

Tip description: straight, very fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	90.0 mm – 3.543 inch
Width	18.0 mm – 0.709 inch
Height	7.0 mm – 0.276 inch

Tip dimensions:

Tip width	0.30 mm – 0.012 inch
Tip thickness	0.30 mm – 0.012 inch