

Product Datasheet

Product code: **372MS.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 372 are Micro scissors with ultra fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, very fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	85.0 mm – 3.346 inch
Width	17.0 mm – 0.669 inch
Height	7.5 mm – 0.295 inch

Tip dimensions:

Tip width	0.05 mm – 0.002 inch
Tip thickness	0.20 mm – 0.008 inch