

Product Datasheet

Product code: **382QMS.9**

Product type: Surgical Scissors

Made in: Germany

Description: The 382Q are Micro scissors with ultra fine blades, ideal for miniature work in confined spaces. The spring-loaded handles provide the scissors with an ergonomic grip, suitable for comfortable execution of multiple cuts. The Inox material is suitable for cleaning. The blades are hardened.

Tip description: straight, very fine

Coating/finishing: Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

| | |
|--------|----------------------|
| OAL | 85.0 mm – 3.346 inch |
| Width | 17.0 mm – 0.669 inch |
| Height | 7.5 mm – 0.295 inch |

Tip dimensions:

| | |
|---------------|----------------------|
| Tip width | 0.15 mm – 0.006 inch |
| Tip thickness | 0.30 mm – 0.012 inch |