

Product Datasheet

Product code: **384**

Product type: UltraCut Scissors

Made in: Italy

Description: The 384 are UltraCut scissors with fine blades, specifically designed for miniature works such as cutting suture fibers. The sharp and spear shape allows to access small areas and guarantees a precise and clean cuts. The micro-serrated blade allows non-slip cutting. The Inox material combined with the polish finishig is the best solution for cleaning. The blades are hardened to 56HRC hardness.

Tip description: straight, fine, sharp, micro-serrated

Coating/finishing: Mirror Polished

Body material: Stainless Steel (Hardened AISI 420)

Picture:



Product dimensions:

OAL	90.0 mm – 3.543 inch
Width	45.0 mm – 1.772 inch
Height	4.0 mm – 0.157 inch

Tip dimensions:

Tip width	0.50 mm – 0.020 inch
Tip thickness	1.00 mm – 0.039 inch