

Product Datasheet

Product code: **394NS**

Product type: UltraCut Scissors

Made in: Italy


Description: The 394 are UltraCut scissors with fine blades, specifically designed for miniature works such as cutting suture fibers. The short cutting edges combined with the extra long handles boost the leverage, making cutting easily and effortlessly. The sharp blades guarantee a precise and clean cuts. The special Inox material combined with the polish finishig is the best solution for cleaning. The blades are hardened to 60HRC hardness.

Tip description: straight, fine, sharp

Coating/finishing: Mirror Polished

Body material: Stainless Steel (Hardened AISI 440)

Picture:



Product dimensions:

OAL	112.0 mm – 4.409 inch
Width	46.0 mm – 1.811 inch
Height	4.5 mm – 0.177 inch

Tip dimensions:

Tip width	0.70 mm – 0.028 inch
Tip thickness	0.70 mm – 0.028 inch