

Product Datasheet


- Product code:** **3CFR.SA.1**
- Product type:** Plastic Replaceable Tip Tweezers
- Made in:** Switzerland
- Description:** Type 3CF Plastic replaceable tip tweezers are suitable to applications which require holding a small component firmly. The Carbon Fiber reinforced ESD Plastic provides the tweezers with ESD and non-scratching properties.
- Tip description:** straight, fine, pointed
- Tip material:** Plastic reinforced with Carbon Fiber (CF)
- Coating/finishing:** Sandblasted
- Body material:** Anti-Magnetic Anti-Acid Stainless Steel (SA)

Picture:



Tweezer dimensions:	OAL	130.0 mm – 5.118 inch
	Width	10.0 mm – 0.394 inch
	Height	16.0 mm – 0.630 inch

Tip details:



Tip dimensions:	Tip width	0.70 mm – 0.028 inch
	Tip thickness	0.70 mm – 0.028 inch