

Product Datasheet

Product code: **3SG.CX.0**

Product type: Biology Tweezers

Made in: Switzerland

Description: Type 3SG high precision tweezers are suitable to delicate applications which require holding a very small component firmly while granting good visibility. The Superalloy material provides improved performance in severe environments. The line serration provides the handles with extra grip, also preventing any possible slips.

Tip description: straight, ultra fine, pointed, superior finish

Coating/finishing: polished

Body material: Anti-Magnetic Anti-Acid Superalloy (CX)


Picture:




Tweezer dimensions:

OAL	120.0 mm – 4.724 inch
Width	10.0 mm – 0.394 inch
Height	12.0 mm – 0.472 inch

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

Tip width	0.08 mm – 0.003 inch
Tip thickness	0.04 mm – 0.002 inch