

## Product Datasheet

**Product code:** **4.DX.0**

**Product type:** Biology Tweezers

**Made in:** Switzerland


**Description:** Type 4 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The high Alloy Austenitic Stainless Steel provides improved performance in severe environments.

**Tip description:** straight, ultra fine, sharp, superior finish

**Coating/finishing:** polished

**Body material:** Anti-Magnetic Anti-Acid High-alloy Stainless Steel (DX)


**Picture:**



**Tweezer dimensions:**

OAL	110.0 mm – 4.331 inch
Width	10.0 mm – 0.394 inch
Height	10.0 mm – 0.394 inch

**Tip details:**



**Tip dimensions:**

Tip width	0.06 mm – 0.002 inch
Tip thickness	0.03 mm – 0.001 inch