

## Product Datasheet

**Product code:** 4.KO.0

**Product type:** Biology Tweezers

**Made in:** Switzerland

**Description:** Type 4 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The Cobalt Alloy has a high fatigue resistance and maintain its elasticity throughout its lifespan, resulting in exceptional durability that is unmatched by other materials.

**Tip description:** straight, ultra fine, sharp, superior finish

**Coating/finishing:** Polished

**Body material:** Cobalt Alloy (KO)


**Picture:**




**Tweezer dimensions:**

OAL	110.0 mm – 4.331 inch
Width	9.5 mm – 0.374 inch
Height	7.0 mm – 0.276 inch

**Tip details:**



Dimensions mm



Dimensions inch



**Tip dimensions:**

Tip width	0.08 mm – 0.003 inch
Tip thickness	0.04 mm – 0.002 inch