

## Product Datasheet

**Product code:** 4SG.CX.0

**Product type:** Biology Tweezers

**Made in:** Switzerland

**Description:** Type 4SG high precision tweezers are suitable to delicate applications which require holding a very small component firmly while granting good visibility. The Superalloy material provides improved performance in severe environments. The line serration provides the handles with extra grip, also preventing any possible slips.

**Tip description:** straight, ultra fine, sharp, superior finish

**Coating/finishing:** polished

**Body material:** Anti-Magnetic Anti-Acid Superalloy (CX)


**Picture:**




**Tweezer dimensions:**

OAL	110.0 mm – 4.331 inch
Width	10.0 mm – 0.394 inch
Height	11.0 mm – 0.433 inch

**Tip details:**



Dimensions mm



Dimensions inch



**Tip dimensions:**

Tip width	0.06 mm – 0.002 inch
Tip thickness	0.03 mm – 0.001 inch