

Product Datasheet

Product code: **5.DX.0**

Product type: Biology Tweezers

Made in: Switzerland

Description: Type 5 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The high Alloy Austenitic Stainless Steel provides improved performance in severe environments.

Tip description: straight, ultra fine, sharp, superior finish

Coating/finishing: polished

Body material: Anti-Magnetic Anti-Acid High-alloy Stainless Steel (DX)


Picture:



Tweezer dimensions:

OAL	110.0 mm – 4.331 inch
Width	10.0 mm – 0.394 inch
Height	11.0 mm – 0.433 inch

Tip details:



Tip dimensions:

Tip width	0.05 mm – 0.002 inch
Tip thickness	0.02 mm – 0.001 inch