

Product Datasheet

Product code: **5.KO.0**

Product type: Biology Tweezers

Made in: Switzerland

Description: Type 5 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting good visibility. The Cobalt Alloy has a high fatigue resistance and maintain its elasticity throughout its lifespan, resulting in exceptional durability that is unmatched by other materials.

Tip description: straight, ultra fine, sharp, superior finish

Coating/finishing: Polished

Body material: Cobalt Alloy (KO)


Picture:




Tweezer dimensions:

OAL	110.0 mm – 4.331 inch
Width	10.0 mm – 0.394 inch
Height	11.0 mm – 0.433 inch

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

Tip width	0.05 mm – 0.002 inch
Tip thickness	0.02 mm – 0.001 inch