

# Product Datasheet

**Product code:** **571.SA.1**

**Product type:** Components Positioning Tweezers

**Made in:** Switzerland

**Description:** Type 572 component positioning tweezers are suitable to delicate applications which require holding a 5-millimeter cylindrical device vertically. The line serration provides the handles with extra grip.

**Tip description:** straight, arc shape, grooved

**Coating/finishing:** Sandblasted

**Body material:** Anti-Magnetic Anti-Acid Stainless Steel (SA)


**Picture:**



**Tweezer dimensions:**

OAL	145.0 mm – 5.709 inch
Width	11.0 mm – 0.433 inch
Height	16.0 mm – 0.630 inch

**Tip details:**



**Tip dimensions:**

Tip width	5.20 mm – 0.205 inch
Tip thickness	1.60 mm – 0.063 inch