

Product Datasheet

Product code: **5B45.DX.0**

Product type: Surgical Tweezers

Made in: Switzerland


Description: Type 5B45 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting maximum visibility. The bent geometry provides maximum visibility and increased ergonomics. The high Alloy Austenitic Stainless Steel provides improved performance in severe environments.

Tip description: bent, very fine, sharp, superior finish

Coating/finishing: Polished

Body material: Anti-Magnetic Anti-Acid High-alloy Stainless Steel (DX)

Picture:



Tweezer dimensions:	OAL	110.0 mm – 4.331 inch
	Width	10.0 mm – 0.394 inch
	Height	11.0 mm – 0.433 inch

Tip dimensions:	Tip width	0.10 mm – 0.004 inch
	Tip thickness	0.06 mm – 0.002 inch