

Product Datasheet

Product code: **5BPSG.S.9**

Product type: Surgical Tweezers

Made in: Germany

Description: Type 5BPSG high precision tweezers are suitable to delicate applications which require holding a component very firmly while granting maximum visibility. The bent geometry provides maximum visibility and increased ergonomics. The line serration provides the handles with extra grip, also preventing any possible slips.

Tip description: bent, very fine, pointed

Coating/finishing: Polished

Body material: Stainless Steel (S)

Picture:



Tweezer dimensions:	OAL	85.0 mm – 3.346 inch
	Width	8.0 mm – 0.315 inch
	Height	10.0 mm – 0.394 inch

Tip dimensions:	Tip width	0.25 mm – 0.010 inch
	Tip thickness	0.20 mm – 0.008 inch