

Product Datasheet

Product code: **7.KO.0**

Product type: Biology Tweezers

Made in: Switzerland


Description: Type 7 high precision tweezers are suitable to delicate applications which require holding a small component firmly while granting maximum visibility. The curved geometry provides maximum visibility and increased ergonomics. The Cobalt Alloy has a high fatigue resistance and maintain its elasticity throughout its lifespan, resulting in exceptional durability that is unmatched by other materials.

Tip description: curved, ultra fine, sharp, superior finish

Coating/finishing: Polished

Body material: Cobalt Alloy (KO)


Picture:




Tweezer dimensions:

OAL	120.0 mm – 4.724 inch
Width	15.0 mm – 0.591 inch
Height	11.0 mm – 0.433 inch

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

Tip width	0.07 mm – 0.003 inch
Tip thickness	0.03 mm – 0.001 inch