

Product Datasheet

Product code: **7SG.CX.0**

Product type: Biology Tweezers

Made in: Switzerland


Description: Type 7SG high precision tweezers are suitable to delicate applications which require holding a very small component firmly while granting maximum visibility. The curved geometry provides maximum visibility and increased ergonomics. The Superalloy material provides improved performance in severe environments. The line serration provides the handles with extra grip, also preventing any possible slips.

Tip description: curved, ultra fine, superior finish

Coating/finishing: polished

Body material: Anti-Magnetic Anti-Acid Superalloy (CX)


Picture:




Tweezer dimensions:

| | |
|--------|-----------------------|
| OAL | 120.0 mm – 4.724 inch |
| Width | 15.0 mm – 0.591 inch |
| Height | 11.0 mm – 0.433 inch |

Tip details:



Dimensions mm



Dimensions inch



Tip dimensions:

| | |
|---------------|----------------------|
| Tip width | 0.06 mm – 0.002 inch |
| Tip thickness | 0.05 mm – 0.002 inch |