

Product Datasheet

Product code: **810.S.P**

Product type: Surgical Tweezers

Made in: PL

Description: Type 810 high precision tweezers are suitable to delicate applications which require holding a component very firmly while granting good visibility. The curve geometry provides good visibility and increased ergonomics. The line serration provides the tips and the handles with extra grip, also preventing any possible slips.

Tip description: curved, fine, beveled, line serrated

Coating/finishing: Sandblasted

Body material: Stainless Steel (S)

Picture:



Tweezer dimensions:	OAL	100.0 mm – 3.937 inch
	Width	7.0 mm – 0.276 inch
	Height	15.0 mm – 0.591 inch

Tip dimensions:	Tip width	0.50 mm – 0.020 inch
	Tip thickness	0.50 mm – 0.020 inch