

Product Datasheet

Product code: 810A.S.P

Product type: Surgical Tweezers

Made in: PL

Description: Type 810 high precision tweezers are suitable to delicate applications which require holding a component very firmly while granting maximum visibility. The bent geometry provides maximum visibility and increased ergonomics. The line serration provides the tips and the handles with extra grip, also preventing any possible slips.

Tip description: curved, fine, beveled, line serrated

Coating/finishing: Sandblasted

Body material: Stainless Steel (S)

Picture:



Tweezer dimensions:

OAL	100.0 mm – 3.937 inch
Width	9.0 mm – 0.354 inch
Height	15.0 mm – 0.591 inch

Tip dimensions:

Tip width	0.50 mm – 0.020 inch
Tip thickness	0.50 mm – 0.020 inch