

Product Datasheet

Product code: **810S.S.P**

Product type: Surgical Tweezers

Made in: PL

Description: Type 810 high precision tweezers are suitable to delicate applications which require holding a component very firmly while granting good visibility. The line serration provides the tips and the handles with extra grip, also preventing any possible slips.

Tip description: straight, fine, beveled, line serrated

Coating/finishing: Sandblasted

Body material: Stainless Steel (S)

Picture:



Tweezer dimensions:

OAL	100.0 mm – 3.937 inch
Width	7.0 mm – 0.276 inch
Height	15.0 mm – 0.591 inch

Tip dimensions:

Tip width	0.50 mm – 0.020 inch
Tip thickness	0.50 mm – 0.020 inch