

Product Datasheet

Product code: **ADSON.SA.B3**

Product type: Surgical Tweezers

Made in: Pakistan

Description: Type ADSON precision tweezers are suitable to applications which require holding a big component very firmly. The line serration provides the tips and the handles with extra grip, also preventing any possible slips.

Tip description: straight, thick, beveled, line serrated

Coating/finishing: Polished

Body material: Anti-Magnetic Anti-Acid Stainless Steel (SA)

Picture:



Tweezer dimensions:	OAL	125.0 mm – 4.921 inch
	Width	14.0 mm – 0.551 inch
	Height	15.0 mm – 0.591 inch

Tip dimensions:	Tip width	1.80 mm – 0.071 inch
	Tip thickness	1.00 mm – 0.039 inch