

Product Datasheet

Product code: **AF4135ZTE**

Product type: Digital Microscopes

Made in: Taiwan

Description: The AF4135ZTE excels with its versatile magnification range of 10x - 330x, ensuring detailed exploration across various applications. This model also includes a polarizer, features FLC, and is compatible with the WF-20 Wi-Fi streamer.

Special features: Polarizer - FLC

Resolution: 1.3 Megapixel (1280x1024)

Magnification: 10-330X

Picture:



Product dimensions:	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch