

## Product Datasheet

**Product code:** **AF4135ZTLE**

**Product type:** Digital Microscopes

**Made in:** Taiwan

**Description:** The AF4135ZTLE stands out with its wireless-ready capability (requires the optional WF-20 adapter), offering 10x - 200x magnification and a 1.3 Megapixel Edge sensor for exceptional clarity. It features an integrated polarizer to reduce reflections, extensive measurement functionality, and Flexible LED Control (FLC) for versatile lighting adjustments.

**Special features:** Polarizer - FLC - LWD

**Resolution:** 1.3 Megapixel (1280x1024)

**Magnification:** 10-200X

**Picture:**



<b>Product dimensions:</b>	OAL	114.0 mm – 4.488 inch
	Width	33.0 mm – 1.299 inch
	Height	33.0 mm – 1.299 inch